



- *Lightweight cannula with or without gas port*
- *Pyramidal tip or conical tip - sharp trocar  
- easy insertion*
- *Silicone valve - airtight seal*
- *Threaded cannula - efficient fascia anchoring*
- *Reduced procedure costs - cost effective*
- *Sterile, single use - no cleaning, no risk of  
cross contamination*



**LiNA Port**  
*Single use 5mm  
Cannula and Trocar System*

**The LiNA Single Use Cannula and Trocar System** meets the demand for quality and efficiency at the most competitive price. The small size allows easy access to the abdomen and good and easy manoeuvrability of instruments.

**2 Lengths - Short and long Ports:** Short 62mm and Long 88mm. Diameter inside cannula: 5,8mm. Allows easy insertion of all 5mm instruments. Suitable for all laparoscopic procedures and all patients.

**Airtight Silicone Double Seal Valve.**  
Secures position of instrument. Easy exchange of instruments. Enables the use of instruments down to 3mm.

**Non-traumatic Fascia Thread.**  
Safe and firm anchoring of cannula.

**Conical Tip or Pyramidal Tip.**  
For easy insertion through fascia.



The LiNA Port is available in two packings:

Single pack



Dual pack with one LiNA Clear Sight Wipe



**Ordering information:**

CE  
LiNA Port 0459

REF.	Product	Packing
1005-1	5mm Cannula and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1005-2	5mm Cannula and Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2005-1	5mm Cannula and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2005-2	5mm Cannula and Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets
1007-1	5mm Cannula and Conical Trocar - short	Single pack, sterile / Box of 5 sets
1007-2	5mm Cannula and Conical Trocar - short	Dual pack*, sterile / Box of 5 sets
2007-1	5mm Cannula and Conical Trocar - long	Single pack, sterile / Box of 5 sets
2007-2	5mm Cannula and Conical Trocar - long	Dual pack*, sterile / Box of 5 sets
1008-1	5mm Cannula with gas port and on/off valve and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1008-2	5mm Cannula with gas port and on/off valve + one Cannula without gas port and Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
1009-1	5mm Cannula with gas port without on/off valve and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1009-2	5mm Cannula with gas port without on/off valve + one Cannula without gas port and Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets

LiNA distributor in your country:

\*2 x Cannula + 1 x Trocar + 1 Clear Sight Wipe

