



Emptying & Drainage Bags



Innovative and practical



SAP by Curas™ minimise risk of cross infection



Safe and convenient disposal of urine



Flexi-Funnel™ ensures excellent closed connection



Suitable for Mitomycin C and cytotoxic waste handling



REF 20460

C3 Emptying Bags

C3 Emptying Bags are specifically designed for emptying urine from Standard Urine Bags, Closed System Bags and Urine Meters. In contrast with any reusable jug or plastic container, it allows the disposal of urine in a hygienic and safe manner. Curas C3 Emptying Bag is highly recommended for demanding and complex patient situations. The C3 bags are genuinely intended to minimise the risk of infection by providing a closed connection during the emptying of the urine.

C2 Short-Term High-Risk Drainage Bags

C2 Standard Bags with a long tube and male connector are developed for safe disposable of urine from, Secretion Bags, Night Bags, Catheters and External Male Catheters.



REF 20467

C4 Closed System Bags (Short-Term High-Risk)

C4 Closed System Bags were developed for the high-risk procedure of bladder instillation and other short term drainages. It is ideal for use in safe handling of high-risk infections and drugs, such as Mitomycin C residue. These bags are particularly well-suited for use in shorter-term operating room-driven cases, where catheter removal typically occurs within the initial 6 -18 hours.

C2 Standard Bags

REF	Description	Capacity (L)	Tube Length (cm)	Connector Type	Non-Return Valve	SAP by Curas™	Inner Box (pcs)	Outer Carton (pcs)
20467	C2 Short-Term High-Risk Drainage Bag, Non-DEHP, PLUS	2.0	110	Male Connector	Yes	Yes	25	200

C3 Emptying Bags

REF	Description	Capacity (L)	Connector Type	Non-Return Valve	SAP by Curas™	Carry Cord/ Cross Hanger™	Inner Box (pcs)	Outer Carton (pcs)
20460	C3 Emptying Bag, Non-DEHP	2.0	Flexi-Funnel™	Yes	Yes	Yes	25	200
20462	C3 Emptying Bag, Non-DEHP	4.0	Flexi-Funnel™	Yes	Yes	Yes	10	40
20465	C3 Emptying Bag, Non-DEHP	1.5	Medium Flexi-Funnel™	Yes	Yes	Yes	25	200
20469	C3 Emptying Bag, Non-DEHP	2.0	Medium Flexi-Funnel™	Yes	Yes	Yes	25	150
20470	C3 Emptying Bag, DEHP-Free	2.0	Flexi-Funnel™	Yes	Yes	Yes	25	150
20471	C3 Emptying Bag, DEHP-Free	2.0	Flexi-Funnel™	Yes	No	Yes	25	300

C4 Closed System Bags (Short-Term High-Risk)

REF	Description	Capacity (L)	Tube Length (cm)	Connector Type (With Sample Port)	Non-Return Valve	SAP by Curas™	Inner Box (pcs)	Outer Carton (pcs)
20490	C4 Closed System Bag (For High Risk Drainage, Bladder Instillation etc.), Non-DEHP	2.0	110	Male Connector	Yes	Yes	10	80
20491	C4 Closed System Bag (For High Risk Drainage, Bladder Instillation etc.), Non-DEHP	4.0	110	Male Connector	Yes	Yes	10	40