

A.4083

A.4084

A.4084NS

sapi med S.p.A.

THE BEAK

THE BEAK DIAGNOSTIC




I "The Beak" sono anoscopi dalla linea assolutamente innovativa nel campo degli strumenti chirurgici per la proctologia. La forma della punta, chiusa e arrotondata, consente l'introduzione nell'ano senza l'ausilio del mandrino. I "The Beak" offrono allo specialista l'opportunità di esporre in modo chiaro e soddisfacente una porzione (un quarto di circonferenza) di canale anale e retto ampollare.

The proctoscopes "The Beak" are instruments with a uniquely innovative design feature. The closed, rounded tip allows for introduction of the scope directly into the anus, needing no separate introducer. The proctoscopes "The Beak" offer the specialist the opportunity of clearly and satisfactorily exposing a portion (one fourth of the circumference) of the anal canal and rectal ampulla.



THE BEAK



REF A.4083

The Beak, anoscopio operativo fenestrato e serigrafato, ideale per tutti i tipi di interventi anorettali. Non contiene lattice. Sterile e monouso.

Sterile: ETO

Confezione: 25 pezzi


REF A.4083

The Beak is a surgical silk-screened proctoscope with a full-length open channel through which anorectal procedures may be performed. Latex free. Disposable&sterile.

Sterilized: ETO

Quantity per box: 25 pieces

THE BEAK DIAGNOSTIC



REF A.4084NS

The Beak Diagnostic, anoscopio fenestrato e serigrafato, ideale per la visita diagnostica ambulatoriale. Non contiene lattice. Monouso.

Confezione: 25 pezzi

REF A.4084NS

The Beak Diagnostic is a silk-screened proctoscope with a full-length open channel. It is ideal instrument for visiting outpatient diagnostic. Latex free. Disposable.

Quantity per box: 25 pieces



REF A.4084

The Beak Diagnostic, sterile e monouso.

Sterile: ETO

Confezione: 25 pezzi

REF A.4084

The Beak Diagnostic, disposable&sterile.

Sterilized: ETO

Quantity per box: 25 pieces

SAPI MED S.p.A.

15121 Alessandria (ITALY)

Tel. +39.0131.348.109 r.a. • Fax +39.0131.348.383

www.sapimed.com • e-mail: info@sapimed.com

Distributed by