


PES Pilot HDpro



STROBOSCOPE



TYMP Module



DIGITAL 5 Module



CHIP-ON-THE-TIP

# Pilot Diagnostic Center

Modular system for diagnostics


**otopront**<sup>®</sup>

# Pilot Diagnostic Center

Modular system for diagnostics

Modularly expandable  
by **PLUG&PLAY**

The **PES Pilot HDpro** can be upgraded to a compact but complete **diagnostic center** right at your workplace. In addition to the **Full HD camera (optionally 4K)**, it can be retrofitted with our proven components such as **LED stroboscopy**, **tympanometry (TYMP)** and **A-mode sonography (DIGITAL 5)**. All findings can be archived in a single patient database.



## PES Pilot HDpro

- + Full HD/4K photo and video recording
- + operation by touch screen or foot pedal
- + patient and analysis software and full connectivity




+



### STROBOSCOPE

- + freely adjustable high-power LED light source
- + Hz and dB data integrated in the video
- + frame-by-frame analysis; image extraction


+



### TYMP Module

- + fully automatic impedance measurement
- + compliance measurement in just two seconds
- + ipsi or contralateral triggering of the stapedius reflex

+



### DIGITAL 5 Module

- + ultrasonic method: A-mode, 4.25 MHz
- + gain infinitely variable (65 - 120 dB)
- + optional single or multiple image display

+



### CHIP-ON-THE-TIP-CAMERA

- + High resolution, bright and true color images
- + FEES and general dysphagia examinations
- + Compatible with disposable covers

## Happersberger otopront GmbH

Medical Equipment Manufacturer  
D-65329 Hohenstein · P.O Box · Germany  
Tel.+49 6120-9217-0 · Fax+49 6120-9217-60  
www.otopront.de · info@otopront.de