

ENDO LINEAR CUTTER STAPLER AND RELOAD UNITS




QUALITY CARE THROUGH INNOVATION



Product Name	Product	Description
Universal P-60 Endo Linear Cutter	CEAG45	Endo Linear Cutter (Standard Handle)
	CEAG60	Endo Linear Cutter (Extra Long Handle)
Articulating Reload Units	CADG-30T	Articulating Reload Unit, 30mm, 2.5mm Staples Height
	CADG-30N	Articulating Reload Unit, 30mm, 3.5mm Staples Height
	CADG-45T	Articulating Reload Unit, 45mm, 2.5mm Staples Height
	CADG-45N	Articulating Reload Unit, 45mm, 3.5mm Staples Height
	CADG-45S	Articulating Reload Unit, 45mm, 4.0mm Staples Height
	CADG-45D	Articulating Reload Unit, 45mm, 4.2mm Staples Height
	CADG-45R	Articulating Reload Unit, 45mm, 4.4mm Staples Height
	CADG-60T	Articulating Reload Unit, 60mm, 2.5mm Staples Height
	CADG-60N	Articulating Reload Unit, 60mm, 3.5mm Staples Height
	CADG-60S	Articulating Reload Unit, 60mm, 4.0mm Staples Height
	CADG-60D	Articulating Reload Unit, 60mm, 4.2mm Staples Height
	CADG-60R	Articulating Reload Unit, 60mm, 4.4mm Staples Height

P-60 ENDO LINEAR CUTTER WITH 60 DEGREES RELOADS

NON-STOP ARTICULATING LEVER



Panther Healthcare Medical Equipment Co., Ltd
Beijing | Suzhou | New York | Sao Paulo | Mumbai | Mexico D.F

www.pantherhealthcare.com
Email: service@pantherhealthcare.com


www.pantherhealthcare.com

P-60 ENDO LINEAR CUTTER WITH 60 DEGREES NON-STOP ARTICULATING LEVER



TECHNOLOGY P-60

- Extra Gripping Technology
- Reduce forward tissue migration
- Lateral diffusion of the fluids
- Improve tissue compression



Automatic pulling
back system



All RUs can pass through
12mm trocar

Non-stop articulation system

Universal handle with better
design and ergonomics for
easy 360° knobbing

±60°

Wide-angle staple
cartridge design
Provides greater operational flexibility

Powerful clamping force
EG design