



New Generation Endo Stapling System

P3H





TECHNOLOGY P₃H

- Lateral diffusion of the fluids
- Stepped cartridge surface
- Better tissue compression

The P₃H technology guarantees better stapling delivering better hemostasis and safety

BENEFITS

- Enhanced anvil delivers better compression & manipulation of the tissue
- Stapling with P₃H technology at the areas of greater tissue tension
- Hemostasis & safety

FEATURES

Stepped
cartridge



Drop-shaped
staple groove



Progressive staple
pushing technology



Enhanced
anvil






● Purple




● Black



● Brown

● Curved Tip Reloads for the most challenging procedures



Stronger clamping force



Four staples at the end of the staple line



Interchangeable cartridge



Coating



New Generation Endo Stapling System

Product Name	Product	Description	Color	Trocar diameter
UNIVERSAL Endo Linear Cutter	CEAB-30	Endo Linear Cutter (Standard Handle)		
	CEAB-45	Endo Linear Cutter (Short Handle)		
	CEAB-60	Endo Linear Cutter (Extra Long Handle)		
Articulating Reloading units	CADD-30EVS	Endo Articulating reloading unit 30 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	 Brown	12mm
	CADD-30ENTS	Endo Articulating reloading unit 30 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple	12mm
	CADD-30EPTS	Endo Articulating reloading unit 30 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black	15mm
	CADD-45EVS	Endo Articulating reloading unit 45 mm for Vascular/ Medium tissue with compression range from 0.88 to 1.5mm	 Brown	12mm
	CADD-45ENTS	Endo Articulating reloading unit 45 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple	12mm
	CADD-45EPTS	Endo Articulating reloading unit 45 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black	15mm
	CADD-60EVS	Endo Articulating reloading unit 60 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	 Brown	12mm
	CADD-60ENTS	Endo Articulating reloading unit 60 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple	12mm
Articulating Reloading Units with Curved Tip	CADD-60EPTS	Endo Articulating reloading unit 60 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black	15mm
	CADE-30EVS	Endo Articulating reloading unit 30 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown	
	CADE-45EVS	Endo Articulating reloading unit 45 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown	
	CADE-45ENTS	Endo Articulating reloading unit 45 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple	