



# 「*Urofino*」 Disposable Electronic Flexible Ureteroscope


Shanghai SeeGen Photoelectric Technology Co., Ltd.

**SEEGEN**  
vision borderless





Less weightness and ergonomic design make the scope procedure easier to control.

Easy to control all scope procedure.



2 light guide lens make the viewing more bright.

Taper-shape distal tip makes insertion easier.



Bending angle 275°

<b>Model</b>	<b>UV-US100-H</b>
<b>Field of view</b>	110°
<b>Depth of view</b>	5-50mm
<b>Deflection up/down</b>	275°/275°
<b>Insertion dia.</b>	3.0mm
<b>Distal end dia.</b>	3.0mm
<b>Working length</b>	750mm
<b>Working channel</b>	1.2mm



Shanghai SeeGen Photoelectric Technology Co.,Ltd.  
 Add.: 3 Floor, Building No.1 4299 Jindu Road, Minhang Dist., Shanghai  
 China Email: endoexp@163.com

CE<sub>0123</sub> EN ISO13485:2016